

RX-ACID SCAVENGER DOSAGE CHART

Compressor Oil Volume

Number of 2 oz. Dosages of Rx-Acid Scavenger
Based Upon Acid Contamination

Gallons (Liters)		Light to Moderate Acid	Severe Acid	After Flushing with Rx11-flush
1/8	(.47)	1 (2 oz. bottle)	1	1
1/4	(.95)	1	1	1
3/8	(1.4)	1	1	1
1/2	(1.9)	1	1	1
3/4	(2.8)	1	2	1
1	(3.8)	1	2	1
1 1/2	(5.7)	2	3	1
2	(7.6)	2	3	1
2 1/2	(9.5)	3	4	2
3	(11.4)	3	5	2
4	(15.1)	4	6	2
5	(18.9)	5	8	3
6	(22.7)	6	9	3
7	(26.5)	7	11	4
8	(30.3)	8	12	4
9	(34.1)	9	14	5
10	(37.8)	10	15	5
11	(41.6)	11	17	6
12	(45.4)	12	18	6

Note: System design as well as refrigerant and refrigeration oil capacities will vary. As a result, it is recommended that the system be tested again after several days of operation to insure that all the acid was scavenged and neutralized. If not, repeat the application.

